

## Aidacare Floorline Bed Accessories



FL250 Bed accessories to improve the functionality of the bed for individual users

- **Self Help Pole:** The Self Help Pole provides manoeuvring and raising support to a resident. The Self Help Pole holder facilitates varied positioning and is securely attached behind the bed head
- **Split folding Side Rails (Pair):** Left/Right pair of moulded rails provides manoeuvring assistance to a resident and/or a safety barrier. Can be positioned at the head end, foot end or a pair at both ends. Folds down to lowest position beside the bed platform
- **Folding Assist Bar:** Provides manoeuvring assistance to a resident. Folds down when not in use, out of the way of transfers
- **IV Pole:** The height adjustable IV Pole can also be manoeuvred within its holder, which is discreetly attached behind the bed head
- **Battery Backup:** Battery Backup is a recommended safety option to prolong bed operation by up to 15 lifts during a power outage
- **Width Extension Kit:** Allows conversion of a Single FL250 into a King Single. Note - King Single FL250 Bed also sold separately
- **EasyMove Transport System:** The EasyMove Transport System is an ideal transport solution designed to minimise manual handling. It allows effortless lifting and tilting of the bed to its transport position, while the large 100 mm castors offer smooth and easy manoeuvring for one person. The EasyMove is also a perfect maintenance platform facilitating easy access to mechanical parts under the bed platform

- ✔ Contracted DVA product
- ✔ Contract approved

## Warranty Information

- 1 Year

| Code      | Display Name   |
|-----------|--|
| BEB046540 | FL250 - Self Help Pole - bed attached (sold sep)                   |
| BEA829140 | FLX / AC4 - Self Help Pole   |
| BEB046525 | FL250 / FLX - Split folding Side Rails (Pair)                      |
| BEB046520 | FL250 - Grab/Assist Rail - Folding Single (each sold sep)          |
| BEA029710 | FLX / AC4 - IV Pole  |
| BEB046530 | FL250 - IV Pole  |
| BEB046590 | FL250 - V1 Battery Backup  |
| BEB046500 | FL250 / FLX - King Single Conversion Kit (8pc set - for retro fit) |
| BEB046235 | FL250 - Acute Conversion Kit                                       |
| BEB046592 | FL250 / FLX - V2 Battery Backup                                    |
| BEB046616 | FL250 / FLX - USB Port Upgrade Kit                                 |
| BEA050055 | FLX - EasyMove Transport System                                    |