

## Elevating Legrest



- A great range of accessories to fit the Aspire Manual Wheelchair range to improve functionality and user experience
- Click [here](#) to ensure this accessory fits your specific wheelchair model

- ✔ Contracted DVA product
- ✔ Contract approved

### Warranty Information

- 1 Year

### Compliance

AS3696-1:2008



| Code      | Display Name                                      |
|-----------|---|
| WCA848100 | Elevating Legrest - Left - 300/320/350/400/410mm  |
| WCA848110 | Elevating Legrest - Right - 300/320/350/400/410mm |
| WCA848120 | Elevating Legrest - Left - 450/460/500/510mm      |
| WCA848130 | Elevating Legrest - Right - 450/460/500/510mm     |
| WCA848140 | Elevating Legrest - Left - 570/610mm              |
| WCA848150 | Elevating Legrest - Right - 570/610mm             |