

## Liko Original Highback Sling




- Can be used to lift most patients
- Good choice for lifting from a horizontal position
- Provides support for the entire back and head
- Enables a semi-reclined sitting posture

### Warranty Information

- 1 Year



 Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS367854	Original High Back MOD 200/210 Net Poly - MEDIUM	200 kg	Net Polyester
LSS367856	Original High Back MOD 200/210 - Net Poly - XLARGE	500 kg	Net Polyester