

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
MWC427046	Wheelchair Rea Azalea Manual Standard Transit 20in (as per script)	\$ 4,219.00	\$ 2,109.50	2
LSS361800	Lifting Chair - Raizer II w/Remote Carry Casem Seat Cover & Battery Charger - 150kg	\$ 8,399.00	\$ 2,519.70	1
CHP197828	Chair Pressure Relief - Monaco Chair - Large - Electric Rise & Recline - 160kg	\$ 5,912.00	\$ 1,773.60	1
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	2
BEB046060	Bed - Aspire Homecare Bed w/Transport Bracket - Single - 4 Function - 180kg	\$ 2,371.00	\$ 1,185.50	2
MWS442600	Wheelchair Folding Push/Pull - 440mm Leverdrive - Left Hand	\$ 4,629.00	\$ 2,314.50	1
MWS442700	Wheelchair Folding Push/Pull - 440mm Leverdrive - Right Hand	\$ 4,629.00	\$ 2,314.50	1
CHP197765	Chair Powerlift Recline - Aspire Signature - Dual Action - Large - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	2
MWC427042	Wheelchair Rea Azalea Manual Standard Transit 18in (as per script)	\$ 4,219.00	\$ 2,109.50	1
MWC427300	Wheelchair Quickie LXI / LX (As per script)	\$ 3,989.00	\$ 1,994.50	1
CHP199700	Chair Powerlift Recline - Air Comfort - Airwing - Sml 500x430x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	1
BEA007215	Alarm - Invisabeam - Remote with Power Supply	\$ 539.00	\$ 161.70	7
LIC350110	Ceiling Hoist - Prism P200 Portable - Manual Traverse w/Quick Release System - 200kg	\$ 3,400.00	\$ 1,020.00	1
LSS390592	Standing & Transport Aid - Shifty Transfer Assist	\$ 3,099.00	\$ 929.70	1
BEB046050	Bed - Aspire COMMUNITY CARE Bed - 5 Function Folding with Tburg - SHPOle - 180kg	\$ 2,877.00	\$ 1,438.50	1
BEB045900	Bed - Aidacare AC3 Single - 3 Function - CL - Beech Head & Foot Boards - 220kg	\$ 2,829.00	\$ 1,414.50	1
PTB533795	WChair Back - JAY J3 - Posterior & Lateral - Shoulder Height Medium 18in	\$ 1,289.00	\$ 644.50	2
BEA007408	Alarm - Cura - Falls Prevention Demo Kit	\$ 2,419.00	\$ 725.70	1
MWC427124	Wheelchair Glide Series 3 (As per script)	\$ 2,369.00	\$ 1,184.50	1
CHP197750	Chair Powerlift Recline - Aspire Signature - Dual Action - Medium - Lagoon - 155kg	\$ 2,269.00	\$ 1,134.50	1

Item Number	Item Description	WAS	NOW	QTY
BEA007215	Alarm - Invisabeam - Remote with Power Supply	\$ 539.00	\$ 161.70	7
BEA007408	Alarm - Cura - Falls Prevention Demo Kit	\$ 2,419.00	\$ 725.70	1
BEA007516	Nurse Call Test Box	\$ 120.00	\$ 36.00	1
BEA007523	Alarm - Cura - Cable Connector - 6.35 Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	10
BEA007539	Alarm - Cura - Pad Connector - USB to 6.35mm R/Angled Plug Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	2
BEA007554	Alarm - Cura - Cordless Multiport Expander (5 devices) - Ring and Sleeve Configured	\$ 287.00	\$ 86.10	1
BEA007850	Alarm - Sofihub TEQ Home (1 x SOFIHUB TEQ Home + 8 Sensors)	\$ 1,419.00	\$ 709.50	1
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	4
BEA024121	Bed End Set - CS7 wide - Drake 42" w/ Ctrl Panel & T-burg - Biltmore Cherry	\$ 789.00	\$ 236.70	1
BEA024126	Bed Extender - Invacare - 100mm - for Invacare All CS Beds - Adjustable Width	\$ 550.00	\$ 165.00	2
BEA028100	Overbed/Chair Table - Aspire - Wide C Base - Thermoform Recessed Split Tilt Top - Cream	\$ 349.00	\$ 174.50	2
BEA029705	Self Help Pole - Aidacare AC/Ward Beds	\$ 184.00	\$ 92.00	1
BEA828150	Bedside Cabinet - 4 Drawer + Hutch (select colour)	\$ 609.00	\$ 304.50	1
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	3
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	2
BEB045900	Bed - Aidacare AC3 Single - 3 Function - CL - Beech Head & Foot Boards - 220kg	\$ 2,829.00	\$ 1,414.50	1
BEB046050	Bed - Aspire COMMUNITY CARE Bed - 5 Function Folding with Tburg - SHPOle - 180kg	\$ 2,877.00	\$ 1,438.50	1
BEB046060	Bed - Aspire Homecare Bed w/Transport Bracket - Single - 4 Function - 180kg	\$ 2,371.00	\$ 1,185.50	2
BEB046490	Bed End Set - Moulded ABS for Acute Care - Aidacare FL250 Bed	\$ 689.00	\$ 206.70	1
BEC046060	Bedside Commode w/ Swing Out Bowl - K Care - 440mm Slate - Phosphate	\$ 382.90	\$ 191.45	2
BEM469737	Mattress - Icare - Latex - Long Single	\$ 1,709.00	\$ 854.50	1
BTC062400	Shower Commode Seat - K Care - H/Made - Closed 500mm	\$ 749.00	\$ 224.70	1
BTC062500	Shower Commode Seat - K Care - H/Made - Closed 550mm	\$ 839.00	\$ 251.70	1
BTC0671105	Shower Commode Elevating Legrest - Aspire - LEFT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC0671155	Shower Commode Elevating Legrest - Aspire - RIGHT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC085100	Bowl Carrier - K Care - 320mm - for K Care Shower Commodes - S/Steel	\$ 354.00	\$ 106.20	2
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front	\$ 461.00	\$ 138.30	2
BTS105400	Bath Board Alton Bariatric 190kg	\$ 285.00	\$ 85.50	1
BTS108305	Bath Transfer Bench - Extension Legs	\$ 67.00	\$ 20.10	1
BTS117700	Shower Stool w/ Arms & Padded Back & Seat - 125kg	\$ 319.00	\$ 95.70	2
BTS123890	Shower Curtain Rings - - Shower curtain rings - Clear	\$ 7.49	\$ 2.25	1
BTS147130	Shower trolley hose kit	\$ 245.00	\$ 73.50	2
BTS151110	Shower Chair - Saljol - Spa - Medium - Dark Grey - 150kg	\$ 220.00	\$ 66.00	2
CHP197055	Chair Powerlift Recline - Pride - C101 Riser Recliner Pecan - 510x510x700mm - 1 Motor	\$ 1,391.50	\$ 695.75	1
CHP197750	Chair Powerlift Recline - Aspire Signature - Dual Action - Medium - Lagoon - 155kg	\$ 2,269.00	\$ 1,134.50	1
CHP197765	Chair Powerlift Recline - Aspire Signature - Dual Action - Large - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	2
CHP197770	Head/Armrest Cover Set - Aspire Signature Recliners - Mink	\$ 46.30	\$ 13.89	1
CHP197772	Head/Armrest Cover Set - Aspire Signature Recliners - Heritage	\$ 46.30	\$ 13.89	2
CHP197828	Chair Pressure Relief - Monaco Chair - Large - Electric Rise & Recline - 160kg	\$ 5,912.00	\$ 1,773.60	1
CHP197835	Armrest Cover Set w/ 1in Padding - Phoenix & Sorrento Chairs	\$ 659.00	\$ 197.70	1
CHP199700	Chair Powerlift Recline - Air Comfort - Airwing - Sml 500x430x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	1
CHP216400	Chair Tray Small - Princess Chair/Bed - Air Comfort	\$ 149.00	\$ 44.70	2
CHP217100	Chair Footrest - Air Comfort - Tilt Bed	\$ 280.90	\$ 84.27	1
CHP217300	Chair Tray Large - Tilt / Deluxe Bed/Daily Chair - Air Comfort	\$ 262.00	\$ 78.60	1
CON244500	Cello Prevail Pants Large (Pkt 18)	\$ 65.00	\$ 19.50	2
CON244550	Cello Prevail Pants Extra Large (Pkt 14)	\$ 65.00	\$ 19.50	4
DIN266600	Tumbler Autoplas 230ml Green	\$ 5.49	\$ 1.65	14
DLG271180	Utensil Holder - Rolyan - Small	\$ 77.00	\$ 23.10	1
DLG271270	Utensil Handles - Foam Tubing - Grey 12mm x 1m (31mm)	\$ 26.00	\$ 7.80	1
DLG272215	Lighting - Outdoor Sensor Light	\$ 71.00	\$ 21.30	1
DLK276110	Cutlery - Kings Double lightweight foam handle w strap	\$ 41.20	\$ 12.36	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
DLK276380	Cutlery - Amefa - Contoured Knife Angled	\$ 83.69	\$ 25.11	1
DLK276390	Cutlery - Amefa - Contoured Fork Angled Left	\$ 52.31	\$ 15.69	1
DLK276400	Cutlery - Amefa - Contoured Fork Angled Right	\$ 52.31	\$ 15.69	1
DLK276410	Cutlery - Amefa - Contoured Spoon Angled Left	\$ 93.00	\$ 27.90	1
DLR279310	EZY Reacher with Safety Lock 20" (S.M.S.)	\$ 35.00	\$ 10.50	1
DLR279320	EZY Reacher with Safety Lock 30" (S.M.S.)	\$ 37.00	\$ 11.10	2
FBS283201LM	Thumb Brace - Rolyan - Collum CMC Thumb Brace - Left - Med	\$ 120.00	\$ 36.00	1
FBS283530	Ankle Brace - BodyAssist - Elastic w/ Heel Loop Anchor - (UNIVERSAL Sz) - Beige	\$ 40.00	\$ 12.00	1
FBS283590	Lumbar - BodyAssist - Thermal Lumbo Sacral Suppt - Blk - Lge-XL	\$ 125.00	\$ 37.50	1
FBS284960	Ankle Brace - McDavid - Thermal Wrap - Blk (UNIVERSAL Sz)	\$ 26.00	\$ 7.80	1
FBS301142	Hand - Glove - Isotoner - Therapeutic pressure gloves - Full Finger - Sml	\$ 63.00	\$ 18.90	1
FBS301144	Hand - Glove - Isotoner - Therapeutic pressure gloves - Full Finger - Lge	\$ 63.00	\$ 18.90	1
FBS301278XXL	Lumbar - Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - XXL	\$ 94.00	\$ 28.20	1
FBS301285XL	Lumbar - Oppo - Light Industrial Back Suppt - XL	\$ 53.40	\$ 16.02	1
FBS301363	Sling - Rolyan - Shoulder Support Abduction Sling - Lge	\$ 161.00	\$ 48.30	1
FBS301379L2	Wrist Brace - Rehband - Manu Comfort T Stable Wrist - Left - Med / - Lge	\$ 120.00	\$ 36.00	1
FBS301421	Wrist Brace - Rolyan - Breathoprene Wrist Splint - Lge - Left	\$ 43.80	\$ 13.14	1
FBS301519L	Knee Brace - Oppo - Neoprene Wrap w/ Patella Tracker Suppt - Lge	\$ 67.00	\$ 20.10	1
FBS301519XL	Knee Brace - Oppo - Neoprene Wrap w/ Patella Tracker Suppt - XL	\$ 67.00	\$ 20.10	1
FBS301526L	Knee Brace - Rehband - Genu Carezza Wrap w/ Stays - O/Patella - Lge	\$ 178.00	\$ 53.40	1
FBS301538M51	Knee Immobiliser - Oppo - Length 51 cm - Med	\$ 65.00	\$ 19.50	1
FBS301541	Knee Immobiliser - Post-op Knee Immobiliser - Med	\$ 50.00	\$ 15.00	1
FBS301569	Ankle Brace - Oppo - Rigid Air Lite w/ straps - Med	\$ 84.00	\$ 25.20	1
FBS301569L	Ankle Brace - Oppo - Rigid Air Lite w/ straps - Lge	\$ 84.00	\$ 25.20	1
FBS303514	Gloves - Smart Glove (each) Large	\$ 76.00	\$ 22.80	2
FPP301900	Reach N Range Pulley with Metal Bracket	\$ 105.00	\$ 31.50	1
GMD300320	Sphygmomanometer - Digital Semi Auto Hand Held	\$ 198.00	\$ 59.40	1
GMD332780	Sphygmomanometer Obese / Thigh Cuff Velcro (S.M.S.)	\$ 32.00	\$ 9.60	2
GMD332955	Sphygmomanometer Cuff & Bladder Set Adult 2 Tube Navy Blue	\$ 32.20	\$ 9.66	1
GMS336050	Stocking Accessories - "Compression Assist" limb lubricant - 60ml bottle	\$ 37.40	\$ 11.22	3
GSS340251	Cushion Cover - Waterproof cover for Latex Multi-use Roll	\$ 58.00	\$ 17.40	1
LIC350110	Ceiling Hoist - Prism P200 Portable - Manual Traverse w/Quick Release System - 200kg	\$ 3,400.00	\$ 1,020.00	1
LIC353380	Ceiling Hoist - Prism - Sonell Charger for Transactive Xtra	\$ 535.00	\$ 160.50	1
LIC355904	Hub Strap Guide - Prism For P200 - Thick White	\$ 26.50	\$ 7.95	2
LSS361600	Sling Quickfit Deluxe Padded Large (S.M.S.)	\$ 361.00	\$ 108.30	1
LSS361800	Lifting Chair - Raizer II w/Remote Carry Casem Seat Cover & Battery Charger - 150kg	\$ 8,399.00	\$ 2,519.70	1
LSS373310	Lifter Sling - Invacare Hammock - Large	\$ 224.00	\$ 112.00	1
LSS374305	Sling - General Purpose - Fabric - XL	\$ 327.00	\$ 163.50	1
LSS376100	Sling - Supreme Hygiene - Fabric Small	\$ 327.00	\$ 163.50	1
LSS377388	Lifter Pivot Frame - 6 point - to suit Oxford Presence Hoist	\$ 288.00	\$ 86.40	1
LSS377395	Lifter Pivot Frame - 4 point - to suit Oxford Advance Hoist	\$ 719.00	\$ 215.70	1
LSS377415	Lifter Yoke to suit Alto Hoist	\$ 205.00	\$ 61.50	7
LSS377746	In-Situ Straight Leg Sling-parasilk GPHS-Loops-Only-Med	\$ 619.00	\$ 309.50	1
LSS377763	Lifter Sling - Allegro - Head Support Pivot Fabric - Large	\$ 444.00	\$ 222.00	2
LSS377768	Lifter Sling - Allegro - Head Support Pivot - Mesh - Medium	\$ 424.00	\$ 212.00	1
LSS377769	Lifter Sling - Allegro - Head Support Pivot - Mesh - Large	\$ 444.00	\$ 222.00	2
LSS383600	Lifter Pivot Frame w/ Keyhole Sling Points - Kerry	\$ 1,189.00	\$ 356.70	1
LSS390592	Standing & Transport Aid - Shifty Transfer Assist	\$ 3,099.00	\$ 929.70	1
LSS390960	Lifter Sling - Aspire - SML - Low Friction General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	1
LSS393540	Spreader Bar - Liko Universal SlingBar 600 - Max 300 kg (660 lbs).	\$ 490.00	\$ 147.00	1
LSS393791	Quick Release Universal Hook & Q-Link Assembly - Liko	\$ 199.00	\$ 59.70	1
LSS394150	Sling - SupportVest w/ padding Mod 91 - Liko - XL - w/2 Loop Clips Fits Sabina II 350 SlingBar.	\$ 749.00	\$ 224.70	1
LSS394160	Sling - SafetyVest Mod 93 - Liko - Med - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 809.00	\$ 242.70	1
LSS394180	Sling - ComfortVest Mod 95 - Liko - Sml - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 779.00	\$ 233.70	1
LSS394190	Sling - ComfortVest Mod 95 - Liko - Med - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 799.00	\$ 239.70	1
LSS394200	Sling - ComfortVest Mod 95 - Liko - Lge - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 829.00	\$ 248.70	1
LSS394210	Sling - ComfortVest Mod 95 - Liko - XL - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 849.00	\$ 254.70	1
LSS394290	Sling (Walking) - MasterVest - Front Close Mod 64 - Liko - Poly - XL - 200kg	\$ 1,019.00	\$ 305.70	1
LSS394800	Liko Extension Loop, green, 12 cm (5 in)	\$ 17.40	\$ 5.22	2
MBT401105	Battery - Topin - 6 Volt 5Amp/Hour - for Weigh Chairs	\$ 31.10	\$ 9.33	2
MWC427042	Wheelchair Rea Azalea Manual Standard Transit 18in (as per script)	\$ 4,219.00	\$ 2,109.50	1
MWC427046	Wheelchair Rea Azalea Manual Standard Transit 20in (as per script)	\$ 4,219.00	\$ 2,109.50	2
MWC427124	Wheelchair Glide Series 3 (As per script)	\$ 2,369.00	\$ 1,184.50	1
MWC427300	Wheelchair Quickie LXI / LX (As per script)	\$ 3,989.00	\$ 1,994.50	1
MWS442370	Wheelchair Folding - Aspire METRO X - Self Prop - METALLIC ORANGE - 100kg ^	\$ 679.00	\$ 339.50	1
MWS442600	Wheelchair Folding Push/Pull - 440mm Leverdrive - Left Hand	\$ 4,629.00	\$ 2,314.50	1
MWS442700	Wheelchair Folding Push/Pull - 440mm Leverdrive - Right Hand	\$ 4,629.00	\$ 2,314.50	1
PED490500	Seating Support - Firefly - GoTo Size 1 w/Advanced Headrest - Blue	\$ 850.00	\$ 595.00	1
PTA493900	HipSaver Slim Fit Veterans Kit Female Medium	\$ 529.00	\$ 370.30	1
PTA494020	HipSaver Slim Fit Veterans Kit Male Large	\$ 529.00	\$ 370.30	1
PTA494153	HipSaver - QuickChange - Pant Only - Female - Med	\$ 65.00	\$ 45.50	2
PTA501310	Heel Holder - #101	\$ 89.00	\$ 62.30	2
PTA508480	WChair Armpad Trough Liner - Pride - Large - Gel	\$ 203.00	\$ 142.10	2
PTA509010	Leg Protectors - Pride - Full Length - Gel	\$ 293.00	\$ 205.10	1
PTA509021	Underarm Crutch Protector - Pride - Gel	\$ 58.00	\$ 40.60	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
PTA509300	Heel Protection - Heelift Suspension Boot Traction	\$ 140.00	\$ 98.00	1
PTA513250	HipSaver - SlimFit - Starter Kit (3-Pants, 1-Pads) - Female - Lge	\$ 269.00	\$ 188.30	2
PTA513800	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Lge	\$ 170.00	\$ 119.00	1
PTA518400	JAY Seat Solid Insert 18x16 Std	\$ 126.00	\$ 88.20	1
PTA518500	JAY Seat Solid Insert 18x18in	\$ 126.00	\$ 88.20	1
PTA518600	JAY Seat Solid Insert 16x16in	\$ 100.00	\$ 70.00	1
PTA520501	Boot - Sheepskin Hi Temp - Tall Slipper - O/Toe - ONE Sz FITS ALL	\$ 276.00	\$ 193.20	1
PTA522200	Footplate Cover Velour (S.M.S.)	\$ 60.00	\$ 42.00	2
PTA534101	Limb/Cast Protection - Limbo - Adult Hand - Waterproof	\$ 51.00	\$ 35.70	1
PTB532200	WChair Back - JAY J2 - Deep Contour 16in	\$ 1,429.00	\$ 714.50	1
PTB532220	WChair Back - JAY J2 - Deep Contour 18in	\$ 1,429.00	\$ 714.50	1
PTB532511	JAY Back Basic 18 x 18in Fixed Clamp	\$ 529.00	\$ 264.50	1
PTB533685	J3 Back Posterior & Deep Contour Shoulder Height Small 18in	\$ 1,529.00	\$ 764.50	1
PTB533717	WChair Back - JAY J3 - Posterior - Mid Thoracic Med 20in	\$ 1,119.00	\$ 559.50	1
PTB533777	WChair Back - JAY J3 - Posterior - Shoulder Height Tall 20in	\$ 1,179.00	\$ 589.50	1
PTB533779	WChair Back - JAY J3 - Posterior & Lateral - Shoulder Height Small 16in	\$ 1,289.00	\$ 644.50	1
PTB533781	WChair Back - JAY J3 - Posterior & Lateral - Shoulder Height Medium 20in	\$ 1,319.00	\$ 659.50	1
PTB533795	WChair Back - JAY J3 - Posterior & Lateral - Shoulder Height Medium 18in	\$ 1,289.00	\$ 644.50	2
PTB533836	WChair Back - JAY J3 - Laterals - Swing Away 6x4in	\$ 759.00	\$ 379.50	1
PTB533886	WChair Back - JAY J3 - Extra Mounting Hardware	\$ 358.00	\$ 179.00	2
PTC544640	Cushion - ROHO Single Compartment HP 23.75x22in (13x12cell)	\$ 1,289.00	\$ 644.50	1
PTC550075	Cushion - Advantage - Otto Bock - 16 x 18in	\$ 904.20	\$ 452.10	1
PTC554500	Cushion - ROHO Harmony 16x16in	\$ 409.00	\$ 204.50	1
PTC561271	Cushion - Jay X2 16 x18in	\$ 759.00	\$ 379.50	2
PTC568200	Cushion - JAY - Care 20x17.5in	\$ 949.00	\$ 474.50	1
PTC571902	Cushion - JAY - Lite 14x14in	\$ 509.00	\$ 254.50	1
PTC573427	Cushion - JAY - Union 22 x 24in	\$ 779.00	\$ 389.50	1
PTC585149	Cushion - JAY - Lite 17x17	\$ 482.00	\$ 241.00	1
PTC591100	Cushion - ROHO Quadtro Select LP 13x15in (7x8cell)	\$ 1,119.00	\$ 559.50	1
PTC594899	Cushion - ROHO Hybrid Elite 1R 15x15in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594956	Cushion - ROHO Hybrid Elite 2R 15x15in	\$ 1,079.00	\$ 539.50	1
PTC596150	Cushion Cover - ROHO Universal HP 18.5x16in	\$ 153.00	\$ 45.90	10
PTC597300	Cushion Cover - ROHO Standard LP 15x16.75in	\$ 153.00	\$ 45.90	1
PTM638015	Mattress - Aspire HYBRID - King Single - 9 Cell - 1980x1050x170mm - 180kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTS640460	Stockings - JOBST UltraSheer - Waist - 15-20mmHg - P/hose - C/Toe - Med - Natural	\$ 124.00	\$ 37.20	1
PTS640920	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Natural	\$ 102.00	\$ 30.60	1
PTS641873SRS	Stockings - GloriaMed - Unisex (Cotton) - Knee - 20-30mmHg - O/Toe - Sml/ reg/short- Bge	\$ 103.00	\$ 30.90	1
RAM65648	Ramp - Scooter Ramp - 1830 x 620 x 710	\$ 1,209.00	\$ 362.70	1
RAM656481	Ramp - Threshold - 4"	\$ 276.00	\$ 82.80	1
RAM656800	Ramp - Aspire Transport - 1520 x 850mm - 400kg	\$ 939.00	\$ 281.70	1
SCA668220	Scooter Acc - Cane Holder Double (Rear) w/flag holder (requires L Piece)	\$ 115.00	\$ 34.50	1
TSE672526-S	Evacuation Sheet - Aidacare Lifecomfort - Single - 1900mm x 720mm - 300kg	\$ 150.00	\$ 75.00	3
TSE673090	Draw Sheet - HD - Blue - 940mm x 156mm	\$ 305.00	\$ 91.50	2
TSE673155	Handi-Lift Walk Belt X Small 58-70cm	\$ 145.00	\$ 43.50	2
WAA686130	Walker Skis White 19mm	\$ 79.00	\$ 23.70	2
WAA686140	Walker Skis White 22mm	\$ 79.00	\$ 23.70	3
WAA686370	Tip Ferrule Bowlers 16mm for Walking Sticks	\$ 36.80	\$ 11.04	1
WAA686375	Tip Ferrule Bowlers 19mm for Walking Sticks	\$ 36.80	\$ 11.04	1
WAA686911	Aspire Walker Accessories - Underseat Basket - XL	\$ 31.80	\$ 9.54	2
WAA686980	Walker Acc - Gutter Walker - Forearm Supports Cuff Pad	\$ 118.00	\$ 35.40	6
WAA695615	Crutch & Cane Holder - Aspire VOGUE INDOOR Walkers	\$ 12.50	\$ 3.75	4
WAF699215	Forearm Walker - Gutter Frame (Oval Tube) - TITAN	\$ 549.00	\$ 164.70	1
WAF704500	Forearm Walker - K Care - Walking Tutor - Child w/ Glides - 600-790mm	\$ 664.40	\$ 199.32	1
WAF705600BL	Seat Walker - Aspire PUSH DOWN - Blue - 7in Wheels - 130kg	\$ 207.00	\$ 103.50	1
WAF709650CH	Seat Walker - Aspire SUPERLITE Walker - Champagne - HBrakes - 6in Wheels - 113kg ^	\$ 219.00	\$ 109.50	2
WAF712100BL	Knee Walker - Aspire - Blue - HBrakes - Bin Wheels - 136kg	\$ 301.00	\$ 150.50	1
WAF751100	Seat Walker - Saljol - Page Indoor Walker - Small - Charcoal - 150kg	\$ 380.00	\$ 114.00	2
WAF751110	Seat Walker - Saljol - Page Indoor Walker - Medium - Charcoal - 150kg	\$ 380.00	\$ 114.00	2
WAS748110	Walking Stick - Freedom - Cobra Cane - 710-940mm - 110kg	\$ 56.00	\$ 16.80	2
WAS748130	Walking Stick - Freedom - Crook Handle - Aluminium w/Black Handle - 770-970mm - 110kg	\$ 26.40	\$ 7.92	1
WAS748140	Walking Stick - Freedom - Crook Handle - Black - 770-970mm - 110kg	\$ 28.50	\$ 8.55	4
WAS748150	Walking Stick - Freedom - T Handle Plastic - Black - 770-970mm - 110kg	\$ 29.90	\$ 8.97	1
WAS748170	Walking Stick - Freedom - Crook Handle - Bronze w/Black Handle - 770-970mm - 110kg	\$ 23.70	\$ 7.11	1
WAS748200	Walking Stick - Freedom - Crook Handle - Silver - 770-970mm - 110kg	\$ 23.70	\$ 7.11	3
WAS748210	Walking Stick - Freedom - Crook Handle - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748240	Walking Stick - Freedom - Folding - Bronze - 800-9000 - 114kg	\$ 28.50	\$ 8.55	3
WAS748250	Walking Stick - Freedom - Offset Handle - Black - 770-970mm - 110kg	\$ 24.00	\$ 7.20	4
WAS748260	Walking Stick - Freedom - Offset Handle - Bronze - 770-970mm - 110kg	\$ 24.00	\$ 7.20	6
WAS748270	Walking Stick - Freedom - Offset Handle - Chrome - 770-970mm - 110kg	\$ 24.00	\$ 7.20	8
WAS748280	Walking Stick - Freedom - Offset Handle - Tartan - 770-970mm - 110kg	\$ 24.00	\$ 7.20	7
WAS748310	Walking Stick - Freedom - Quad Cane - Large Flat Base - 770-980mm - 110kg	\$ 61.00	\$ 18.30	1
WAS748410	Walking Stick - Freedom - Swan Neck - Leaves - 770-970mm - 110kg	\$ 28.50	\$ 8.55	4
WAS748420	Walking Stick - Freedom - Swan Neck - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAS748430	Walking Stick - Freedom - T Handle Ergo - Black - 770-970mm - 110kg	\$ 24.00	\$ 7.20	5

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
WAS748450	Walking Stick - Freedom - T Handle Ergo - Chrome - 770-970mm - 110kg	\$ 25.20	\$ 7.56	6
WAS748460	Walking Stick - Freedom - T Handle Ergo - Leaves - 770-970mm - 110kg	\$ 24.00	\$ 7.20	6
WAS748470	Walking Stick - Freedom - T Handle Ergo - Tartan - 770-970mm - 110kg	\$ 24.00	\$ 7.20	4
WAS748500	Walking Stick - Freedom - T Handle Ergo Plastic - Chrome - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748520	Walking Stick - Freedom - T Handle Ergo Plastic - Leaves - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAS748550	Walking Stick - Freedom - T Handle Plastic - Chrome - 770-970mm - 110kg	\$ 28.50	\$ 8.55	3
WAS748580	Walking Stick - Freedom - T Handle Plastic - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAT760800	Wash Bowl - Green Autoclavable - 240mm x 110mm - 2000ml	\$ 21.40	\$ 6.42	1
WAT806300	Instrument Trolley - Kerry - 600x500 - 1 Drawer S/S - Epoxy PC Frame	\$ 829.00	\$ 414.50	1
WAT815200	FlexiCart Options - Drawer Dividers with Identification Labels Set 2	\$ 283.00	\$ 84.90	1
WAT816650	FlexiCart Options - Folding Shelf Heavy Duty with Flat Top	\$ 719.00	\$ 215.70	1
WAT818500	FlexiCart Options - Waste Bin (No Bracket)	\$ 105.00	\$ 31.50	1
WAT820180	IV Pole - Aspire - 5 Star Polyurethane Base	\$ 188.00	\$ 94.00	2
WCA848160	WChair Elevating Legrest - Advanced - Left - 320-410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS&RX	\$ 183.00	\$ 54.90	1
WCA848170	WChair Elevating Legrest - Advanced - Right - 320-410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS&RX	\$ 183.00	\$ 54.90	1
WCA848200	WChair Elevating Legrest - Articulating - Left - 570/610mm - Aspire Assist2/Evoke2	\$ 183.00	\$ 54.90	1
WCA848210	WChair Elevating Legrest - Articulating - Right - 570/610mm - Aspire Assist2/Evoke2	\$ 183.00	\$ 54.90	1
WCA848530	WChair Depth Extension Kit - 50mm - Black - Aspire Rehab RS/Rehab RX	\$ 25.80	\$ 7.74	1
WCA848540	WChair Depth Extension Kit - 100mm - Black - Aspire Rehab RS/Rehab RX	\$ 32.30	\$ 9.69	1
WCA848550	WChair Therapist Kit for Hand Adjustment - Aspire Rehab RS/Rehab RX	\$ 38.30	\$ 11.49	1
WCA858008	WChair Stump Support - Right - Heavy Duty WChair (20-24in)	\$ 340.00	\$ 102.00	1
WCA858221	Basket Bracket for Jet 3 Ultra	\$ 70.00	\$ 21.00	1
WCA858390	WChair Elevating Legrest - Aspire EVOKE 400/450/500/HD - Black Right	\$ 120.00	\$ 36.00	1
WCA858460	WChair Anti Tip Bar - Aspire EVOKE	\$ 52.00	\$ 15.60	1
WCA858590	WChair Belt w/ Buckle - Attachment	\$ 93.00	\$ 27.90	2
WCA858822	Headrest Clamp Standard (New Style) - Black	\$ 112.00	\$ 33.60	4
WCA858823	Headrest Clamp Heavy Duty (New Style) - Black	\$ 268.00	\$ 80.40	1
WCA859000	Channel Forearm Pad Small - Ottobock	\$ 239.00	\$ 71.70	1
WCA859170	Modular Hand Pad Horn Small Left	\$ 125.00	\$ 37.50	1
WCA859190	Modular Hand Pad Horn Small Right	\$ 125.00	\$ 37.50	1
WCA859340	Clinical Mounting Kit 52-62cm Black	\$ 356.00	\$ 106.80	1
WCA859370	Clinical Mounting Kit 34-40cm Black	\$ 356.00	\$ 106.80	1
WCA859425	Start Stump support Left	\$ 423.00	\$ 126.90	1
WCA859445	Start Stump support Right	\$ 423.00	\$ 126.90	2
WCA859450	Start Height Adjustable Side Panel - Silver	\$ 323.00	\$ 96.90	1
WCA859496	Start Padded Hip Belt - Adult - 2" x 48" w/ auto style buckle	\$ 97.00	\$ 29.10	1
WCA870236	Power Wheelchair Accessory - Whill C1 - Arm Covers - Red - Right	\$ 218.00	\$ 65.40	1
WCA870300	Power Wheelchair Accessory - Whill C1 - Storage Basket	\$ 246.00	\$ 73.80	1